

DENT-JECT®

Dental Injection Needles

- ▶ Both British and Metric thread are available
- ▶ Triple sharpened cannula tip with advanced siliconized coating for smooth penetration
- ▶ Thin wall design is available for increasing the flow
- ▶ Hub color code according to ISO 6009



Increased flow
Enlarged lumen

Code	Dimensions	Color Code
DN3012	0.30x12mm 30Gx1/2"	Yellow
DN3021	0.30x21mm 30Gx5/6"	Yellow
DN3025	0.30x25mm 30Gx1"	Yellow
DN3035	0.30x35mm 30Gx1 2/5"	Yellow
DN3042	0.30x42mm 30Gx1 2/3"	Yellow
DN2712	0.40x12mm 27Gx1/2"	Grey
DN2716	0.40x16mm 27Gx5/8"	Grey
DN2721	0.40x21mm 27Gx5/6"	Grey
DN2730	0.40x30mm 27Gx1 1/4"	Grey
DN2735	0.40x35mm 27Gx1 2/5"	Grey
DN2740	0.40x40mm 27Gx1 1/2"	Grey

Dental Injection Needle Luer-Lock

Sizes
27G/0.4mm
30G/0.3mm

Packaging
100pcs per box



Dental Paste Needles

- ▶ Blunt round tip for comfortable usage in the gingival sulcus
- ▶ Bendable without significant decrease of the flow rate

Code	Dimensions	Color Code
DPN2012	0.90x12mm 20Gx1/2"	Yellow
DPN2021	0.90x21mm 20Gx5/6"	Yellow
DPN2025	0.90x25mm 20Gx1"	Yellow
DPN2112	0.80x12mm 21Gx1/2"	Green
DPN2121	0.80x21mm 21Gx5/6"	Green
DPN2125	0.80x25mm 21Gx1"	Green
DPN2312	0.60x12mm 23Gx1/2"	Blue
DPN2321	0.60x21mm 23Gx5/6"	Blue
DPN2325	0.60x25mm 23Gx1"	Blue

