

# BERPU-CARE<sup>®</sup>

## Safety Insulin Pen Needles

- ▶ Own patent products
- ▶ Automatically trigger the safety function after injection
- ▶ Achieve safety function without more operation than normal pen needles
- ▶ Self-made sharpened needle and top lubrication technology giving minimum penetration force
- ▶ Uniform hub screw thread can fit with a wide range of insulin pens in the market



Code	Dimensions	Color Code
SIN3005	0.30x5mm 30Gx5mm	
SIN3006	0.30x6mm 30Gx6mm	
SIN3008	0.30x8mm 30Gx8mm	
SIN2905	0.33x5mm 29Gx5mm	
SIN2906	0.33x6mm 29Gx6mm	
SIN2908	0.33x8mm 29Gx8mm	

