


BERPU-CARE®

Safety Scalp Vein Set

- ▶ Exquisite sharp needle for easy puncture and less pain
- ▶ Soft wings for the convenience of fixation
- ▶ Colorful wings for needle size identification
- ▶ Safety shield can be activated to cover the needle immediately against accidental needle stick injury

Code	Dimensions	Color Code
SV2720	0.40x20mm 27Gx3/4"	
SV2620	0.45x20mm 26Gx3/4"	
SV2520	0.50x20mm 25Gx3/4"	
SV2420	0.55x20mm 24Gx3/4"	
SV2320	0.60x20mm 23Gx3/4"	
SV2220	0.70x20mm 22Gx3/4"	
SV2120	0.80x20mm 21Gx3/4"	
SV2020	0.90x20mm 20Gx3/4"	
SV1920	1.10x20mm 19Gx3/4"	
SV1820	1.20x20mm 18Gx3/4"	



Scalp Vein Set

- ▶ Self-made exquisite designed needle, sharp and easy puncture
- ▶ Soft wings for the convenience of fixation
- ▶ Color-coded wings for easy gauge identification
- ▶ Tube length can be customized made

Code	Dimensions	Color Code
V2720	0.40x20mm 27Gx3/4"	
V2620	0.45x20mm 26Gx3/4"	
V2520	0.50x20mm 25Gx3/4"	
V2420	0.55x20mm 24Gx3/4"	
V2320	0.60x20mm 23Gx3/4"	
V2220	0.70x20mm 22Gx3/4"	
V2120	0.80x20mm 21Gx3/4"	
V2020	0.90x20mm 20Gx3/4"	
V1920	1.10x20mm 19Gx3/4"	
V1820	1.20x20mm 18Gx3/4"	

