


# BERPU-JECT®

## Blood Collection Systems

- ▶ Series of blood collection combinations for selection
- ▶ High quality rubber sleeve has excellent self-sealing performance
- ▶ PVC tube length can be made according to customer's requirement
- ▶ Available package: Single blister, single pouch or single PE bag



HLSV2120



HLLSN2138



LLSV2120



LSSN2138



HLS



HLL



HLF



LLV2320



HLLV2320

Component Description	Adapter Luer Slip	Adapter Luer Lock	Holder with Luer Slip Adapter	Holder with Luer Lock Adapter	Holder with Female Adapter	SVS	Safety SVS	Hypodermic Needle	Safety Hypodermic Needle	Outer Diameter (G)	Needle Length (mm)
Product Code	LS	LL	HLS	HLL	HLF	V	SV	N	SN	19-27	20-38
HLS			●								
HLL				●							
HLF					●						
HLSV2120			●			●				21	20
HLSSV2120			●				●			21	20
HLSN2138			●					●		21	38
HLSSN2138			●						●	21	38
HLLV2120				●		●				21	20
HLLSV2120				●			●			21	20
HLLN2138				●				●		21	38
HLLSN2138				●					●	21	38
LSV2120	●					●				21	20
LSSV2120	●						●			21	20
LSN2138	●							●		21	38
LSSN2138	●								●	21	38
LLV2120		●				●				21	20
LLSV2120		●					●			21	20
LLN2138		●						●		21	38
LLSN2138		●							●	21	38